

## INTRODUCTION

- Drug trial for NASH cirrhosis necessitates patients to be classified as F4 by NASH-CRN staging.
- However, F4 classification lacks granularity and a detailed assessment of zonal fibrosis.
- To address this limitation, the approach of Single Harmonic Generation/Two Photon Excitation (SHG/TPE) microscopy-based qFibrosis assessment (qF), enables a fully quantitative evaluation.
- The primary goal of this study is to use qF to compare distribution of zonal fibrosis between two distinct patient groups: those with lean body mass (BMI<25) and those who are overweight/obese (BMI≥25).
- The study aims to enhance our understanding of zonal fibrosis patterns in patients with varying BMI statuses, using qF based analysis.

## METHODS

- n=133 patients were evaluated from Phase 2b Belapectin NASH cirrhosis drug trials (NCT04365868).
- To assess the liver biopsies, a fully quantitative fibrosis assessment method based on NASH-CRN parameters, known as qF, was employed.
- These patients were divided into distinct groups: the Baseline (BL) category included both lean individuals (n=9) and overweight/obese individuals (n=124), while the End-of-Treatment (EOT) category consisted of lean individuals (n=7) and overweight/obese individuals (n=126).
- The statistical analysis of the data was performed utilizing the Wilcoxon rank sum test.
- To visually represent the findings, heat maps were employed, with any statistically significant differences (p≤0.05) indicated in red.

## RESULTS


- At BL, the group exhibited statistically significant periportal fibrosis differences between lean and overweight/obese patients.
- Similar patterns persisted in EOT group, suggesting distinct periportal fibrosis distribution linked to BMI (Figure 1).
- Although the lean group had fewer patients compared to overweight/obese, the consistent observation of significant fibrosis distribution difference in periportal zones in both BL and EOT groups underscores its significance and reproducibility.
- These observations warrant reporting and suggest promising directions for additional exploration, based on BMI categories.

## DISCUSSION


- The ability of SHG/TPE based qF analysis to provide consistent differentiation of periportal fibrosis parameters between lean and overweight/obese patients demonstrates its advantage when compared to conventional NASH-CRN staging.
- As seen in Figure 2, there is more periportal fibrosis in the lean patient biopsy (blue circle) as compared to the obese/overweight patient biopsy (white circle) in both the BL and EOT categories.
- Also in Figure 2, in pericentral zone, although there is similar fibrosis in the lean (blue circle) and obese/overweight (white circle) biopsy in BL, the EOT images depict lesser pericentral fibrosis in the obese/overweight patient (white circle) biopsy compared to lean patient (blue circle) biopsy.
- These observations also indicate histopathological variations in zonal fibrosis distribution in NASH cirrhosis, linked to BMI status.
- Larger patient cohorts should be engaged for deeper insights and validation of these initial findings.

## CONCLUSION

- The potential utility of the outcomes of this study lies in designing NASH drug trials with tailored fibrosis resolution endpoints, accounting for patients' BMI variations.



**Figure 1.** Collagen parameters differ between lean and overweight/obese patients. SHG/TPE microscopy assesses 28 collagen parameters per zone, correlating BL and EOT cohorts via heat map (p≤0.05 in red).



**Figure 2.** SHG/TPE images in lean versus overweight/obese patients in the periportal and pericentral zones.

## ACKNOWLEDGEMENTS

All authors participated in the development of this poster and approved the final version of this poster.

## CONTACT INFORMATION

Dean Tai, PhD  
dean.tai@histoindex.com