

POSTER WED-374

Utilizing AI digital pathology for quantitative measurement, grading and differentiation of lobular and portal inflammation in MASH: preliminary analysis in liver histology

EASL CONGRESS
7-10 May 2025
Amsterdam, the Netherlands



Kutbuddin Akbary¹, Ren Yayun¹, Carolin Lackner², Daniela S. Allende³, Feng Liu⁴, David E. Kleiner⁵, Dean Tai¹

Contact information

Dean Tai, PhD
dean.tai@histoindex.com

[1] HistoIndex Pte Ltd, Singapore; [2] Institute of Pathology, Medical University of Graz, Graz, Austria; [3] Pathology Division, Cleveland Clinic, Cleveland, OH, USA; [4] Peking University Hepatology Institute, Peking University People's Hospital, Beijing, China; [5] Laboratory of Pathology, Centre for Cancer Research, NCI, MD, USA

Introduction

- Histological grading of inflammation is essential for evaluating treatment response in MASH clinical trials
- Differentiating lobular from portal inflammation by light microscopy can be subjective and variable
- AI digital pathology has the potential to offer objective, reproducible assessment through image analysis
- qFibrosis, qSteatosis, and qBallooning were validated using Second Harmonic Generation/Two Photon Excitation Fluorescence (SHG/TPEF) microscopy. However, qInflammation (qi) requires development on H&E images, as SHG/TPEF microscopy has limitations in visualizing inflammatory cells, and qi is more complex to develop and validate

Aim

- To develop and validate an AI-based qiInflammation algorithm using H&E-stained liver biopsies
- To quantify and distinguish lobular and portal inflammation using objective morphometric features using the qi-algorithm

Method

- **Cohorts Utilised:** 269 liver biopsies from different sources (MASH clinical trials and clinical study samples)
 - Randomly divided into training (n=180) and validation (n=89) cohorts
 - Inflammatory foci were classified as described by Brunt et al [1] into portal and intra-acinar (lobular) inflammation ranging from grades 0-3 were used to define Inflammation zones and grades
- **Image Processing:**
 - H&E-stained slides scanned using Leica Aperio scanner
 - Digital images from Haematoxylin and eosin imaging channels of the scanner were processed to segment nuclei, detect inflammatory cells, and identify inflammatory foci
- **Portal Region Detection and Inflammatory Foci Classification:**
 - Portal tracts localized via convolutional neural networks deep learning model trained on annotated blocks
 - Inflammatory foci were classified as lobular or portal based on spatial overlay with portal regions
 - Morphometric parameters, such as density, perimeter, and eccentricity were used to build qi-lobular and qi-portal indices using multiple linear regression on training cohorts and validated on the validation cohort
- **Statistical Analysis:**
 - Spearman correlation with pathologist grades
 - Concordance with pathologist based annotation on the validation cohort measured by Sensitivity (percentage of pathologist-identified inflammation also identified by qi) and Positive Predictive Value (PPV; percentage of qi-identified inflammation also identified by pathologist)
 - Concordance with pathologist inflammation grades also calculated using weighted Kappa, Spearman's correlation, and the area under the receiver operating characteristic (AUROC) analysis.

Results

- **Pathologist Grading Distribution in n=269 biopsies:**
 - Portal inflammation: I0 (14%), I1 (48%), I2 (36%), I3 (2%)
 - Lobular inflammation: I0 (2%), I1 (44%), I2 (47%), I3 (7%)
- **qi detection Performance:**
 - qi algorithm successfully distinguishes and quantifies portal and lobular inflammation from H&E images (Figure 1)
 - Lobular inflammation: Sensitivity 83%, PPV 58%
 - Portal inflammation: Sensitivity 89%, PPV 72%
- **Concordance with Pathologist Grading (Figure 2 and Table 1):**
 - qi-lobular index: Spearman's correlation = 0.556 (p<0.001), Weighted Kappa = 0.32, AUROCs = 0.80-0.88
 - qi-portal index: Spearman's correlation = 0.606 (p<0.001), Weighted Kappa = 0.38, AUROCs = 0.77-0.94




Figure 1: Biopsy image without annotations (Panel A) and the same image annotated using qiInflammation algorithm (green annotations, Panel B) were both graded as **Grade 2 Portal Inflammation** by the pathologist and qiInflammation. Biopsy image without annotations (Panel C) and the same image annotated using qiInflammation algorithm (dark blue annotations, Panel D) were both graded as **Grade 3 Lobular Inflammation** by the pathologist and qiInflammation




Figure 2: Validation cohort data for (A) qi-lobular and (B) qi-portal indices and their correlations with pathologist-based grading of inflammation

Table 1: Performances of the qi-lobular and qi-portal indices. Cutoff values were determined by Youden index. PPV: positive predictive value; NPV: negative predictive value

| A qi-lobular index | | | | | | |
|--------------------|-------|---------|-------------|-------------|------|-----|
| qi-lobular index | AUROC | Cut-off | Sensitivity | Specificity | PPV | NPV |
| I0 vs I1/2/3 | 0.88 | 1.37 | 82% | 100% | 100% | 11% |
| I0/1 vs I2/3 | 0.80 | 1.55 | 67% | 73% | 74% | 64% |
| I0/1/2 vs I3 | 0.82 | 1.83 | 67% | 86% | 20% | 96% |

| B qi-Portal index | | | | | | |
|-------------------|-------|---------|-------------|-------------|-----|------|
| qi-Portal index | AUROC | Cut-off | Sensitivity | Specificity | PPV | NPV |
| I0 vs I1/2/3 | 0.77 | 0.99 | 75% | 46% | 90% | 27% |
| I0 vs I1/2/3 | 0.85 | 1.30 | 74% | 78% | 68% | 83% |
| I0 vs I1/2/3 | 0.94 | 1.69 | 100% | 84% | 13% | 100% |

Conclusions

- Preliminary results show good correlation with pathologist grading and acceptable inter-rater agreement
- Ongoing improvements will focus on refining nuclear segmentation, reducing false positives, and identifying immune cell subtypes
- qi training and validation data completes the qFIBS panel (qFibrosis, qiInflammation, qBallooning, qSteatosis), supporting comprehensive digital assessment of MASH histology for clinical trials

References

1. Brunt et al. Am J Gastroenterol. 1999;94(9):2467-74



Scan to download the poster



EASL CONGRESS
7-10 May 2025
Amsterdam, the Netherlands

DOI: 10.2255/ps0.eu.EASL2025.2025

Metabolism, alcohol and toxicity
Kutbuddin Akbary

WED-374
EASL2025